


Operational immediately


The C 48 series low speed 4 x 400 ml centrifuges offer an excellent capacity/space ratio and enhanced ergonomics: motorized lid opening, 10 program memories for direct access, AWELight illuminated handle to indicate end of run...



C 48
Ventilated



C 48-R
Refrigerated



C 48 / -R	PERFORMANCES
Max. Capacity	4 x 400 ml
Max. Speed	4500 rpm
Max. RCF	3820 xg
Capacity	76 x 5/7 ml blood collection tubes
	28 x 15 ml conical tubes
	12 x 50 ml conical tubes

Classic centrifuges C 48/-R

Specifications	C 48	C 48-R
Max. speed Swing-out rotor Angle rotor	4500 rpm 4400 rpm	4500 rpm 4400 rpm
Max. RCF Swing-out rotor Angle rotor	3820 xg 3030 xg	3820 xg 3030 xg
Max. capacity Swing-out rotor Angle rotor	4 x 400 ml 32 x 15 ml	4 x 400 ml 32 x 15 ml
Control system	Microprocessor	
Speed range	500 to 4500 rpm	
Speed precision	20 rpm	
Temperature Range	-	-9°C to 40°C
Timer	0 - 99h50 and continuous	
Acceleration / Braking Rates	1 / 1	
Programs	10	
Lid Opening	Motorized	
Motor type	Induction	
Certification	EC marked : IEC 1010-1 ; IEC 1010-2-020 Directives	
Dimensions (HxWxD)	280 x 420 x 520 mm	280 x 420 x 750 mm
Crating (HxWxD)	460 x 610 x 650 mm	460 x 610 x 890 mm
Weight (Net/Crated)	50 / 56 kg	70 / 79 kg

Cat. numbers	C 48	C 48-R
230V - 50Hz version	320002001	320003001
120V - 60Hz version	320002011	320003011

Cat. numbers	Accessories
330002001	SO-250 swing-out rotor 4 x 400 ml
330002003	N-250 Round Bucket for SO-250 Swing-out rotor - Set of 4
330002004	BL-250 Lid with bio-containment for N-250 Bucket - Set of 4
330002005	Insert for 14 x 5-7 ml collection tubes - Set of 4
330002009	Insert for 15 x BD 5-7 ml collection tubes - Set of 4
330002010	Insert for 19 x Terumo 5-7 ml collection tubes - Set of 4
330002006	Insert for 12 x 10 ml collection tubes - Set of 4
330002007	Insert for 7 x 15 ml conical tubes - Set of 4
330002008	Insert for 3 x 50 ml conical tubes - Set of 4
330000019	A 48-32.15 Fixed angle rotor 32 x 15 ml tubes



YOUR AWEL PARTNER