

Heating Mantle



Available in various sizes to fit round-bottom flask from 300 to 2000 ml. The snug-fit design of the heating surface, which is in close contact with the glass flask, ensures good thermal efficiency. Max. temperature up to 400°C with accurate temperature control $\pm 1^{\circ}\text{C}$

Model	Capacity	Watts @ 220-240 V
HMIC-F30	300 ml	160 W
HMIC-F50	500 ml	200 W
HMIC-F100	1000 ml	330 W
HMIC-F200	2000 ml	450 W