

Water Bath

Characteristic: static coated, full stainless, stable & precision



Main Data:

Type	DZKW-D-1 single	DZKW-D-2 double	DZKW-D-4 four/line	DZKW-D-6 six/line	DZKW-S-4 four /two lines	DZKW-S-6 six/two lines	DZKW-S-8 eight/two lines
Power (w)	300	500	1000	1500	1000	1500	2000
Volt (v)	220v50Hz	220v50Hz	220v50Hz	220v50Hz	220v50Hz	220v50Hz	220v50Hz
T (°C)	≤±1	≤±1	≤±1	≤±1	≤±1	≤±1	≤±1
Range (°C)	≤±1	≤±1	≤±1	≤±1	≤±1	≤±1	≤±1
Control (°C)	RT+10~ 100	RT+10~ 100	RT+10~ 100	RT+10~ 100	RT+10~ 100	RT+10~ 100	RT+10~ 100
Sensitive (°C)	≤±1	≤±1	≤±1	≤±1	≤±1	≤±1	≤±1
Tolerance (°C)	≤±2	≤±2	≤±2	≤±2	≤±2	≤±2	≤±2
Size (mm)	160×170×90	325×170×90	650×170×90	940×170×90	325×330×90	480×330×90	650×330×90